

Data Sheet of Urine Bag(Leg Bag)

Urine Bag(Leg Bag)	
Material	Non-woven fabric, PVC, etc.
Color	Grey, white
Sterilization	Sterile or non-sterile
Brand	Greaturo



PRODUCTION DESCRIPTION

Greaturo Urine Bag (Leg Bag) is 1500ml size, comes with a needleless sample port and optional valve. It is designed to collect the urine draining from the catheter. This product is used for draining accumulated urine and is especially suitable for the disabled, paralyzed and bedridden patients. This product, compared to regular urine bags, features strap holes and non-woven fabric material, thus serving both as a urine bag and a leg bag.

SPECIAL FEATURES

- With needleless sample port for more safety.
- Available in folding valve, T-valve or single reversal valve.
- With strap holes, dual use for urine bag and leg bag.
- Large 1500ml capacity to collect more urine.

SPECIFICATION

REF	Specification
GCU202003	1500ml, 8.8mm, with needleless sample port, with folding valve, sterile
GCU202004	1500ml, 8.8mm, with needleless sample port, with single reversal valve, sterile
GCU202005	1500ml, 8.8mm, with needleless sample port, with T-valve, sterile
Remark: Non-sterile is available.	



PACKAGING DETAILS

Paper-Plastic Pack	1pc
Master Carton	100pcs

CLASSIFICATION: Class I, sterile or non-sterile

DIRECTIONS FOR USE

1.1 How to Use

- Always wash your Hands before and after emptying the bag.
- Hang bag at or below the patient's bladder level.
- Close outlet before connecting to patient. Replacing the urine bag is at the discretion of the medical specialist staff.
- Periodic observation of this system should be made to ensure urine is flowing freely.
- After collection is full, open outlet valve and drain the contents into collection bottle.
- Then close and clean the outlet valve. Use a disinfection wipe for cleaning and remove any residue at the drain valve.

1.2 Cautions

- For single use only. Dispose of in the waste after use. The product must not be reused. Reuse brings the risk of infection and may impair the performance of the product, up to malfunction.
- This product should only be used under the supervision of medical staff who have received technical training.
- Do not use, if packaging is opened or damaged, if the product is damaged or the expiry date is exceeded.
- During use, the bag must hang upright; there must be no loops and kinks in the tube.
- Do not store in direct sunlight, at extreme temperatures or at high humidity.

1.3 Expiration Date

- The shelf life of the Urine Bag(Leg Bag) is 5 years from the date of manufacture.



(All product, product specifications, data and related descriptions are subject to change without notice to comply with regulatory, improve reliability, function or design or otherwise.)

Rev. UB001 03/24

