


Irrigation Equipments [Evacuation by Irrigation]

■ Irrigation Set ■ Stopper for Irrigation ■ Sleeve for Irrigation

ALCARE considers safety and ease of use in the process of designing products and provides a total series of appliances for comfortable colon irrigation.

■ Irrigation Set



COLOCLEAN PC

Irrigation Set

- This set contains all the appliances needed during and after colon irrigation.
- All parts of the appliances can be put into the attached carrying bag, which is convenient for traveling.

Type	Product Code No.	Specification (Contents of Set)	Contents Per Box
Set	13401	IRRIGATION BAG : 1 pce 2-STEP CONE : 1 pce FACE PLATE : 1 pce DRAINABLE SLEEVE : 5 pcs Carrying bag : 1 pce MINI POUCH-K SUPER : 2 pcs YOU CARE C : 2 pcs	1 set

[Parts]

Product Name	Product Code No.	Specification	Contents Per Box
IRRIGATION BAG	13411	Capacity : 2000ml (with a roll clamp)	1 pce
2-STEP CONE	13412	Silicone (with connector)	1 pce

Product Name	Product Code No.	Specification	Contents Per Box
FACE PLATE	13413	Max Available Size 55mm	1 pce
DRAINABLE SLEEVE	13414	Bag Size (LXW) : 900mm	1 pce

■ Stopper for Irrigation