



PRODUCT CATALOGUE





About Navabharat

Founded in 1934 as a broadsheet Hindi newspaper in Nagpur, NavaBharat Group has earned generations and decades of trust. Currently published from over 14 cities and with sister newspapers in Marathi and English we have one of the strongest readership bases in Central India. While constantly working for the socio economic development of the region, NavaBharat group diversified into real estate development, paper manufacturing & healthcare.

We believe that health is the key to prosperity in every human life. And care is the backbone of health, so we exist to put care in the healthcare.

Nubeno was inceptioned to focus on transformative innovation to deliver smarter, more personalized care for patients and providers. We exist to span the spectrum of healthcare, which reflects in the solutions we create and cater, ranging from Infusion therapy, Respiratory care. We provide answers with precision, delivering breakthrough technologies to prevent, diagnose and treat health needs. Adapting quickly to changes in the world around us, harnessing leading-edge science and technology to deliver the best possible solutions for some of the world's most important health challenges forms our everyday routine.

The entire development process – from the initial idea to being used in the patient's home or as part of the daily hospital routine – we identify possible risks of therapies in dialog with medical specialists and patients. We evaluate them and develop prevention strategies in order to increase treatment quality and safety. With structured risk management, we reduce risks and improve safety in the healthcare field.

Based at the Zero mile of the country, an area of 2 lac Sq. ft. is felicitated with world class facilities of moulding, clean room, sterilization lab, Quality lab under one roof to help enhance outcomes, lower health care delivery costs, increase efficiencies, improve health care safety and expand access to health.



MISSION

Our mission is to improve the quality of life and the quality of our work with our great products and solutions.



VISION

Our vision is to provide solutions and products that are effective. Our vision pushes us to develop high quality solutions.

Infusion Therapy

Infusion therapy involves the administration of medication through a needle or catheter. It is prescribed when a patient's condition is so severe that it cannot be treated effectively by oral medications, the medications are introduced to the vein through a catheter. Diseases commonly requiring infusion therapy include infections that are unresponsive to oral antibiotics, cancer and cancer-related pain, dehydration, gastrointestinal diseases or disorders which prevent normal functioning of the gastrointestinal system, and more.

At Nubeno,, we provide a complete array of the products that aid in Infusion Therapy. These products help in making IV therapy more convenient and safer for the healthcare providers and less painful for the patients. Our entire product range is CE certified and complies with International Quality Standards.



NUFIX IV CANNULA

STERILE EO



Product Features

- Bevel type needle for easier less traumatic insert.
- FEP catheter with X-Ray opaque line.
- Special Flexible and kink resistant catheter is double tapered With rounded tip for smooth introduction and minimum Discomfort to the patient.
- Injection port with non return silicon valve and color coded cap For intermittent and safe administration of drugs During infusion.
- Flexible wings provide easy and secure fixation and also Prevent skin contact with luer lock connector.
- High transparent flash back chamber for an easy view Of returning blood.
- Blood stopper
- Luer lock cap for closing the proximal end of the catheter After removing the needle.

Option Available:

- With Injection Port & Wings
- FEP Clear Catheter
- FEP Radio opaque Catheter
- PTFE Catheter & PU Catheter

Technical Specification

Size	Catheter Length ID mm OD mm		Catheter Length (mm)	Water Flow Rate (ml/min)	Qty. Per Box (pcs.)	Qty. Per Carton (pcs.)
14G Orange	1.7	2.1	45	270	50/100	1000
16G Grey	1.3	1.7	45	180	50/100	1000
17G White	1.1	1.5	45	125	50/100	1000
18G Green	0.9	1.3	45	80	50/100	1000
20G Pink	0.8	1.1	32	54	50/100	1000
22G Blue	0.6	0.9	25	33	50/100	1000
24G Yellow	0.5	0.7	19	20	50/100	1000
26G Violet	0.4	0.6	19	15	50/100	1000



NUBENO THREE WAY STOP COCK

Common Features

- Crystal clear transparent channel
- Minimum priming volume required for accurate drug administration
- Rotating male luer and two female luer connectors
- Arrow indication on the handle to indicate the direction of flow
- Connection meets 6% taper requirement to ensure leakage free fitment
- For pressure up to 4.5 bar (65 psi)

STERILE EO



Packaging

Customization available



Hypodermic Syringe

The word 'Hypodermic' is a medical term that refers to anything below the skin. Case in point, Hypodermic Syringes are syringes with a hollow needle that pierces the skin and injects medicaments into the bloodstream. Likewise, these syringes can also be used to extract blood from the human body. We at Nubeno, offer multiple varieties of Hypodermic Syringes.

The Single-use Hypodermic Syringes are one of the most in-demand single-use syringes in the Indian market.



NUSAFE SYRINGE

Luer slip tip

Our luer slip tip provides a friction-fit connection that requires you to push and twist the syringe tip into the needle hub.

Product Features

- High clarity transparent barrel molded by using premium quality medical grade polypropylene.
- Thermo Plastic Elastomer (TPE) Gasket for smooth movement of the plunger.
- Needles are manufactured by using imported Cannula with precise tip sharpness and strong bonding strength with the hub.
- Needles have an ultra-thin wall permitting rapid and efficient drug delivery.

Technical Specification

Syringe size	1ml	2ml	3ml	5ml	10ml	20ml	50ml
Needle Size 1 " or 1.5" Inch	24G	23G	22G	24G	21G	21G	W/O Needle

Packaging

Syringes Size	Needle Size	Inner Box	Blister Pack Outer Box	Ribbon Pack Outer Box
1ml/cc	24G-28G	100 Pcs	2400 Pcs	2000 Pcs
2ml/cc	21G-26G	100 Pcs	2400 Pcs	2000 Pcs
3ml/cc	21G-26G	100 Pcs	2400 Pcs	2000 Pcs
5ml/cc	21G-24G	100 Pcs	1800 Pcs	2000 Pcs
10ml/cc	21G-24G	100 Pcs	1200 Pcs	1000 Pcs
20ml/cc	21G-22G	50 Pcs	600 Pcs	500 Pcs

STERILE EO



Hypodermic Needle

A Hypodermic (hypo – under, dermic – the skin) needle commonly used with a syringe to inject medicaments into the body or extract fluids from it. We at Nubeno manufacture a variety of high precision Hypodermic Needles designed to offer maximum protection and uplifting results during medical operations. These needles are manufactured keeping in mind the mandates of global industry standards.



Product Features

- Extra thin wall
- Cannula made from imported cold-rolled stainless steel strips in-house
- The tip of needles are grounded using high-precision CNC machines by skilled operators
- Siliconized to minimise drag during penetration and withdrawal
- Vision Camera check for burs, hook and blockage; ensures 100% automated quality control enabling superior sharpness to minimise pain
- Assembled and packaged on totally automated machines

STERILE EO



NUBENO HYPODERMIC NEEDLE

Technical Specification

Colour Code of Hub	Gauge	OD (MM) ± 0.05MM
White	16G	1.6MM
Pink	18G	1.2MM
Cream	19G	1.1MM
Yellow	20G	0.9MM
Green	21G	0.8MM
Black	22G	0.7MM
Blue	23G	0.6MM
Medium Purple	24G	0.5MM
Orange	25G	0.5MM
Brown	26G	0.4MM

Packaging

Blister Packed



NUBENO IV INFUSION SET (VENTED/NON VENTED)

Product Features

- Sharp piercing spike for easy insertion in I.V container.
- Manufactured from non-toxic and transparent material.
- Infusion set with built in air vent and bacteria barrier hydrophobic filter.
- Specially designed transparent cylindrical collapsible drip chamber to visualize the flow rate.
- Disc type 15µm fluid filter at the bottom of drip chamber filters any particulate matter in the I.V fluid.
- 150 cm long super smooth kink resistant tubing.
- Efficient roller controller ensures better flow control.
- Provided with super smooth kink resistant tube for unobstructed flow.
- Fixed Luer lock for secure connection to all standard devices.

STERILE EO



Packaging

- individually packed
- 600 Pcs / Carton

Blood Management System

Blood collection or blood sampling is a crucial aspect in diagnosing various diseases. It not only helps in finding our information on the lack or excess of various elements in the body but also it is the initial step in diagnosing any disease in the body. The medical facilities and hospitals have various choice of blood collection systems. However, there are three major ways through which blood collection is done from the patients, amongst which, Venepunctures is one of the most routinely performed procedures. Evacuated tubes, Multi sample needles, Needle Holder are some of the equipments used in blood collection system.

We at Nubeno offer the range of products which are required in Pre-Analytical Phase of Blood Collection. The products are most important and regularly used equipments in any medical facilities.



NUBENO BLOOD TRANSFUSION SET

Product Features

- For transfusion of blood or blood components
- Transparent & flexible drip chamber
- Approximately 20 drops/ml
- Soft and kink resistant PVC tubing With 200 micron blood filter
- Light weight ball facilitates Auto shut mechanism
- Available with priming filter
- DEHP Free

Option Available :

Vented / Non- Vented / With "Y" Connect

Packaging

Customization available



NUBENO MV SET 100/110/ 150ML

Product Features

- Cylinder, flexible, clear and calibrated measured volume chamber in size 100/110/150ml.
- Chamber has an in built Injection port and air inlet with bacteria filter.
- Burette chamber is made of bio-compatible medical grade transparent styrene polymer, suitable for infusing all types of fluids.
- Unique floating auto shut off valve acts as a floating indicator and automatically shuts off the drain path to prevent the passage of air in the line.
- 15 micron fluid filter in the drip chamber.
- Equipped with a micro drip with reduced drop size of 60 drops/ml. Soft and kink resistant PVC tube, 150 cm long.

Option Available:

- Choice of Latex flush bulb
- "Y" injection port
- Latex free.

Packaging

- Individually packed
- 120 Pcs / Carton

STERILE EO





NUBENO LUNG EXERCISER

Product Features

The Device "Lung Exerciser" used by adults and children to strengthen their respiratory muscles and increase lung volume. The Breather improves the degenerative affects of neuromuscular disease on the respiratory system. It strengthens the muscles associated with breathing and can improve the cough effort.

- Disposable Single patient use.
- Non-Sterile.
- Helps to develop, improve and maintain respiratory fitness.
- Useful in restoring and maintaining lung capacity of post operative ambulatory and COPD patients.
- Air filter for enhanced protection.
- Transparent, well calibrated chambers for accurate estimation of exercise capacity.

Packaging

- Single Unit Pack duplex box.
- Outer Carton- 30 pcs

STERILE EO



Nubeno is one new name in providers and manufacturers of healthcare solutions today. Every service we provide incorporates the entirety of our knowledge and skills, the company's deep understanding of users' needs, and showcase extensive expertise. With its constantly growing portfolio of effective medical care solutions, Nubeno makes a substantial contribution towards protecting and improving people's health.

Our products and services are essential building blocks of healthcare

At Nubeno, our mission is to save and sustain lives. Our Product range under Infusion Therapy, Anaesthesia, Urology, Gastroenterology, Gynaecology, Respiratory Care, Blood Management systems, and more are focused on transformative innovation to deliver smarter, more personalized care for patients and providers. Every day our mission to save and sustain lives is at the forefront as we strive to meet our goal of top-quartile performance. Our objective is to accelerate growth by expanding Nubeno's impact across the continuum of care supported by the evolution of our new product lines. Our company's efforts to achieve top quartile results relative to industry peers and other comparators across three dimensions.

- Patient safety and Quality
- Growth through innovation
- Industry-leading performance

We are dedicated to ensuring that Nubeno is there when and where patients need care, from hospitals and clinics to homes in rural areas and major cities Through:



PRODUCT INNOVATION

Improve sustainability and performance of products and services

Nubeno develops products that save and sustain lives worldwide, delivering on our mission and representing our most important contribution to society. Our commitment to patient safety and quality is integral to innovations that improve the experiences of patients and caregivers who use our products. By integrating sustainability into our product development approach, we consider the full life cycle of our products from design through end-of-life.



PATIENT SAFETY AND QUALITY

Patient safety and product quality are fundamental to fulfilling our mission and ongoing success. We have extensive policies and processes to ensure excellence in these areas, and we continue to innovate new products and services that enhance patient safety and drive improved patient outcomes. We take a holistic, innovative approach to quality across the product life cycle. Building on years of insight, our approach to patient safety and product performance is driven by our commitment to continuous improvement.

INNOVATION THAT EXPANDS ACCESS TO CARE

Guided by our mission and inspired by patients and caregivers who rely on our products, we are dedicated to nurturing, developing and delivering innovative healthcare solutions. Nubeno invest in new collaborations and partnerships with world-renowned institutions, deploying our people, products and grants globally to increase access to Nubeno's lifesaving products.

INNOVATING FOR SUSTAINABILITY ACROSS THE PRODUCT LIFE CYCLE

We are beginning to transition from a linear to a circular approach that maximizes the value gained from materials, while improving environmental performance and the customer experience.

Sustainable Design

We use Environmental, Health, Safety and Sustainability reviews and life cycle assessments to understand product environmental performance and requirements and drive ongoing improvements across the value chain.

Materials Use

We innovate to reduce materials use in our products, decrease waste in manufacturing and packaging, and minimize use of materials of concern.

Patient Safety and Quality

We have an unrelenting dedication to patient safety and quality that is foundational to the company's mission and ongoing success.

Product End-of-Life

We work to repair and reuse electromechanical medical products when possible and collaborate with partners around the world to recycle medical waste.

Product Use

We design products to reduce energy demands and waste during use by customers, caregivers and patients.



NUBENO HEALTHCARE PVT. LTD.

📍 Navabharat Bhawan,
Chhatrapati Square, Wardha Road,
Nagpur - 440015. (M.S.) India

🏭 10 A, Sector No. - 23, Mihan,
Nagpur - 441108. (M.S) India.

☎ 0712 - 2284001 ○ www.nubeno.in

✉ info@nubeno.in