



OPTICLUX 10-2 C T1 (WOOD-LIGHT) D15816000



Illuminance	6 000 lux at 0.15 m
Color Rendering Index CRI (Ra)	95
Color Rendering Index R9	91
Color Rendering Index R13	98
Color temperature	6 500 K
Lens diameter	160 mm
Lens diopter	3.5 dpt (approx. 1.9x magnification)
Lens material	Polymethyl methacrylate (PMMA)
Number of LEDs	105
Operating device	Electronic converter
Power cord	Direct power supply
Power consumption	12 W
Rated input voltage	100 - 240 V
Frequency range	50 / 60 Hz
Luminaire housing	Polyamide (PA)
Weight (Net)	9.1 kg
Operation	Membrane keypad on luminaire head
Dimming	Stepless
Joint mechanism	Spring (counterbalance)
Protection class	I
Degree of protection	IP 20
Standards	EN 60601-1
Special features	continuous rotation of the head (360° swivel axis), scratch-resistant plastic lens, switchable to UV light (Wood Light, 365 nm)
Mounting	Ceiling mount, Ceiling tube can be cut on site (length: 40-116 cm)
Accessories	Intermediate ceiling adapter
Order No.	D15816000-00752633

Light Distribution



Technical Drawing

