


Patient Interface



The key to improving adherence to treatment.

Patient Interface

Our mask solutions for at home and in hospital.



Requirements are highly varied for ventilation masks used in the Homecare area (sleep medicine and outpatient ventilation). The masks have to offer wearing comfort and a perfect fit. They should also be small, lightweight and easy to use. The exhalation system should be quiet with a diffuse flow and a ball-and-socket joint should provide freedom of movement. Finally, the mask should look good too.

The demands on masks used in the hospital have a different focus. There the ventilation masks need a good seal and a stable fit to be able to withstand high pressures. That's particularly important because the pressure difference between inspiration and expiration is often considerable. Depending on their specific uses, the masks have to be reusable (suitable for hygienic reprocessing) or disposable after one use (Single-Patient Use).


Choosing the right mask is crucial for ensuring patient adherence to treatment and its success. That is why our extensive range of masks offers the right model for every requirement.

The product range in the homecare segment includes the JULIA family (without a forehead support) and the internationally award-winning CARA family.

In the hospital segment, the Single-Patient-Use mask JOYCEclinic Full Face is available, while the LENA ventilation mask can be hygienically reprocessed and reused in all settings.

Thanks to this variety, every patient can find exactly the right mask.

The right mask for every requirement:



An overview of our masks



with
double lip



CARA

Nasal mask

vented



CARA
Full Face

Full face
mask

vented



JULIA

Nasal mask

vented



JULIA
Full Face

Full face
mask

vented



JOYCEone

Nasal mask

vented



JOYCEone
Full Face

Full face
mask

vented

non-
vented

Size/ Order number	XS	25590	—	—	—	—	—	—
	S	25600	25630	26710	26740	25260	25290	25270
	M		25640	26720	26750			
	L	25610	25650	26730	26760			
	XL	25580	25680	—	—			
		25590	—	—	—	—	—	—
Pressure range		4 – 30 hPa	4 – 25 hPa	4 - 30 hPa	4 - 35 hPa	4 – 30 hPa	4 – 30 hPa	4 – 30 hPa
Noise level (db(A) at 10 hPa) Sound pressure level		18	19	18	19	19	19	—
Sustainable		—	—	●	●	—	—	—
Ball-and-socket joint		●	●	Flexible mask tube	●	●	●	—
Can be reprocessed for change of patient		●	●	●	●	●	●	●
Made in Germany		●	●	●	●	●	●	●




JOYCEeasy		JOYCEeasy Full Face		JOYCEeasy next Full Face		LENA		JOYCEclinic Full Face	
Nasal mask		Full face mask		Full face mask		Full face mask		Full face mask	
vented	vented	vented	non-vented	vented	non-vented	non-vented	non-vented + AAV	vented	
—	—	—	—	—	—	—	—	—	—
25510	25910	25109	25159	26460	26960	26805 / 26857*	26845 / 26869*	26840 / 26866*	
25520	25920	25119	25169	26470	26970	26815 / 26858*	26855 / 26870*	26850 / 26867*	
25530	25930	25129	25179	26480	26980	26825 / 26859*	26865 / 26871*	26860 / 26868*	
—	—	—	—	—	—	—	—	—	—
4 – 30 hPa	4 – 30 hPa	4 – 40 hPa	4 – 40 hPa	4 – 35 hPa	4 – 35 hPa	4 – 50 hPa	4 – 50 hPa	4 – 50 hPa	
17	13	13	—	12	—	—	—	—	—
—	—	—	—	—	—	—	—	—	—
•	•	•	—	•	—	—	—	—	—
•	•	•	•	•	•	—	—	—	—
•	•	•	•	•	•	•	•	•	•

with
double lip

* in a set of 10

The fit & sizes



Accessories



NV Elbow

In addition to the standard vented variant, the Elbow Full Face NV Set is also available for most full face masks.



NV+AAV Elbow

For the JOYCEclinic Full Face, we offer a non-vented version with an anti-asphyxia valve (AAV), as well as standard vented and non-vented elbows.



Endoscopy Adapter

Endoscopic examination during NIV with our solution for professional users. On many of our masks, simply replace the mask elbow with the Endoscopy Adapter.



Quick-release Cord

The quick-release cord is used in place of headgear clips. Patients with arthritis or limited mobility (e.g., after a stroke) benefit from the simple use.



XL Headgear

This headgear is the top choice of patients with a head circumference of more than 60 cm. It can be combined with all standard masks (except for JULIA).



WILAsilent 2

This external passive exhalation valve is used in clinical environments with NV masks.

	NV Elbow	NV+ AAV Elbow	Endoscopy Adapter NV	Quick-release Cord	XL Headgear (without clips)
JULIA	—	—	—	—	—
JULIA Full Face	—	—	—	—	—
CARA	—	—	—	—	●
CARA Full Face	—	—	—	25623	●
JOYCEone	—	—	—	—	●
JOYCEone Full Face	25275	—	15864	—	●
JOYCEeasy	—	—	—	—	●
JOYCEeasy Full Face	25275	—	15864	—	●
JOYCEeasy next Full Face	25275	—	15864	—	●
LENA	15970	—	15968	26464	●
JOYCEclinic Full Face	25269	26834	15864	Included in standard supply schedule	—



LÖWENSTEIN

medical

Sales + Service

Löwenstein Medical
Arzbacher Straße 80
56130 Bad Ems, Germany
T. +49 2603 9600-0
F. +49 2603 9600-50
info@loewensteinmedical.com
loewensteinmedical.com



Löwenstein Medical Technology
Kronsaalsweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
info@loewensteinmedical.de
loewensteinmedical.com



p10106en2605

© Copyright protected.
Copying or reproduction of any sort requires the explicit approval of Löwenstein Medical.